## Cu(II)-Imprinted Styrene-Divinylbenzene Beads as a New Sorbent For Flow Injection-Flame Atomic Absorption Determination of Copper

Anna Tobiasz, Stanisław Walas\*, Bartosz Trzewik, Przemysław Grzybek, Małgorzata M. Zaitz, Marta Gawin, Halina Mrowiec

Jagiellonian University, Faculty of Chemistry, ul. Ingardena 3, 30-060 Kraków, Poland

## Abstract:

The paper reports preparation and analytical features of a new Cu(II)-imprinted polymer, based on salen-OMe ligand 2,2'-[ethane-1,2diylbis(nitrilo(*E*)methylylidene)]bis(6-allyl-4-methoxy-phenol) and styrenedivinylbenzene matrix, as well as its application to on-line preconcentration and flame atomic absorption determination of copper. Sorbent beads (average diameter of 60-80 µm) were obtained using suspension polymerization technique and employed as a column filling. Copper sorption was the most effective at pH 6.8, whereas the highest elution effectiveness was observed when 0.5 % HNO<sub>3</sub> was applied. The sorbent exhibited good long-term stability and acid resistance. Enrichment factor (EF) of 12 was found for 60 s loading time and loading flow rate of  $4 \text{ mL} \cdot \text{min}^{-1}$ . EF value may be further increased by expanding the loading time and/or flow rate. Batch sorbent capacity in optimal pH conditions was found to be  $0.16 \text{ mmol} \cdot \text{g}^{-1}$  (9.55 mg  $\cdot \text{g}^{-1}$ ) of a dry

polymer. Calcium(II) turned out to be the only significant interferent. Cadmium(II), silver(I), nickel(II), zinc(II) in concentrations lower than about  $1 \text{ mg} \cdot \text{L}^{-1}$  did not disturb preconcentration. Different copper(II) calibration methods such as: set of standards method (SSM), standard addition method (SAM) and combinatory calibration method (CCM) were employed for copper(II) determination in tap water, spring mineral water and certified reference material. Analysis of EU-H-3 reference material confirmed good accuracy of the proposed method. Relative standard deviation (RSD) was 3.2 for standard addition method and 2.8 % for set of standard calibration method. Detection limit for sample consumption 16 mL was 1.03 and 1.07 µg/L respectively.

## Key words:

copper determination, preconcentration, ion imprinted polymer, flow injection, flame atomic absorption spectrometry, salen